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
**(GB) Guedel airway Set**

REV 2022-01

Ningbo Huakun Medical Equipment Co., Ltd.  
No. 208 Binjiang Road  
Cixi High-tech Industrial Development Zone  
Cixi City  
315301 Zhejiang, P.R. China

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## **GB** Contents

|    |  |      |    |
|----|--|------|----|
| 1  | General .....                            | Page | 9  |
| 2  | Safety and Responsibilities .....        | Page | 9  |
| 3  | Modification of the product .....        | Page | 10 |
| 4  | Before using the product .....           | Page | 10 |
| 5  | Using the product .....                  | Page | 11 |
| 6  | Product Maintenance / Reprocessing ..... | Page | 11 |
| 7  | Accessories .....                        | Page | 11 |
| 8  | Replacement parts and repairs .....      | Page | 11 |
| 9  | Disposal .....                           | Page | 11 |
| 10 | Usage .....                              | Page | 12 |
| 11 | Explanation of symbols .....             | Page | 13 |

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## 1 General

- This product may only be used by trained professionals under medical instruction.
- The product is used to ensure breathing in unconscious patients by keeping the upper airways open and preventing the tongue from sliding back into the throat.  
The product also serves as a bite block and bite protection.
- The product is mainly used in emergency care.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

## 2 Safety and Responsibilities

**IMPORTANT! READ CAREFULLY BEFORE USE**  
**KEEP FOR FUTURE REFERENCE**

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- There should be sufficient experience in the use of the product in order to avoid injury to patients.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- The product draws air from the immediate surrounding and must therefore not be used in a toxic atmosphere.
- Use of the product in the presence of oxygen or nitrous oxide, during procedures involving lasers or electrosurgical probes, may cause the product to ignite. This can result in burns and poisoning of the patient.
- Do not apply excessive force at any time during use to avoid injury.
- A fully functioning spare of the product should always be kept ready to hand where the product is being used.
- Do not use silicone-based lubricants.

- The product must not come into contact with pointed or sharp objects to avoid cracks or other damage.
- Select the correct airway size.  
If it is too small, it can cause displacement of the upper airways.  
If it is too large, the epiglottis can be pushed in front of the laryngeal inlet, thus closing the laryngeal inlet.
- Due to the risk of aspiration, the tube must be removed immediately if the patient begins to gag or cough.
- **Single use products must not be reused!**  
**They must be disposed of immediately after use in accordance with chapter 9.**  
The reuse of single use products represents a potential risk for the patient or the user. It can result in contamination and/or reduced functionality.
- Contamination and/or limited functionality of the product can lead to injury, illness or death of the patient.

### **3 Modification of the product**

- Apart from the described maintenance instructions, the product must not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

### **4 Before using the product**

- Do not expose the product to direct sunlight or other light sources. They must be stored and transported in a dry place and in their original packaging.
- Select the correct airway size for the patient.  
The airway size is determined by the distance between the earlobe and the corner of the mouth (see figure a).
- Remove all delivered parts from the packaging.  
If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 10.
- In case of doubt perform a comparison check with a brand-new product.

→ **Damaged products must not be used!**

They must be disposed of in accordance with the instructions in chapter 9.

- The product contains PP.
- Please note that the product is supplied in a non-sterile condition.
- The hygienic preparation of the product lies within the clinical responsibility of the user! The plastic product can be disinfected with a suitable cleaning agent (e.g. alcohol) before use. Sterilization is not permitted!
- The product is now ready for use.

## **5 Using the product**

- Proceed according to the currently valid rules for emergency ventilation.
- Rotate the Guedel airway along its longitudinal axis so that the lower opening of the product faces the palate and insert it into the patient's opened mouth (see figure **b**).
- Insert the product approximately halfway into the oral cavity and rotate it 180° around the longitudinal axis. Continue to insert it slowly so that it follows the anatomical airway.
- Make sure that only the shield of the airway protrudes from the mouth (see figure **c**).
- Make sure the tip of the tongue is positioned against the teeth and is not pushed back into the throat.

## **6 Product maintenance / Reprocessing**

- Single use products must not be reused!
- After termination of the application, the product must be disposed of according to chapter 9.

## **7 Accessories**

- The user is responsible for the correct combination with other medical products.

## **8 Replacement parts and repairs**

- In case of functional or material problems, contact customer service immediately.


## **9 Disposal**

- The product must be disposed of after use.
- Please note the expiry date if applicable.  
The product must not be used after the expiry date!
- Safe disposal is carried out in accordance with current national regulations.


D

10 Verwendung / Stosowanie / Usage


a)


b)

c)

## 11 Symbolerklärung / Objasnienie symboli / Explanation of symbols

|   |  |  |
|--|---|---|
|  Deutschland Kennzeichen   | Symbol Niemiec  | Symbol for Germany  |
|  Polen Kennzeichen   | Symbol Polski   | Symbol for Poland   |
|  Großbritannien Kennzeichen  | Symbol Wielkiej Brytanii  | Symbol for Great Britain  |
|  Gebrauchsanweisung beachten                                       | Przestrzegać instrukcji użytkowania   | Observe the operating instructions  |
|  Warnhinweis   | Wskazówka ostrzegawcza  | Warning   |
|  Hersteller  | Producent   | Manufacturer  |
|  EU-Repräsentant   | Przedstawiciel UE   | EU Representative   |
|  CE-Symbol   | Symbol CE   | CE symbol   |
|  Nicht verwenden wenn Packung beschädigt ist                       | Nie użytkować, jeśli opakowanie jest uszkodzone                                   | Do not use if package is damaged  |
|  Zur einmaligen Verwendung / nur zur Verwendung an einem Patienten | Do jednorazowego użytku / tylko do użytku na jednym pacjencie                     | Single-use only / Only to be used on one patient                                  |
|  Medizinprodukt  | Wyrob medyczny  | Medical Device  |
|  Importeur  | Importer  | Importer  |