Wolfram Droh GmbH, Dammweg 11, 55130 Mainz, Germany

WDI med Sp. z o.o.

ul. Pastelowa 3

70-704 Szczecin

POLSKA







Tel. +48 91 463 05 71 • Fax +48 91 463 04 78

info@wdimed.pl • www.wdimed.pl

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(GB) Splint

REV GA / 2022-07

Anping Guardian Medical Equipment Co., Ltd.
500 meters, Nan Sucun, Anping County
Hengshui City, Hebei Province, 053600
P. R. China

MedNet EC-REP GmbH
Borkstrasse 10
48163 Muenster
Germany

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1 General

- The product may only be used by trained and qualified personnel.
- The product is intended to be used for immobilization and splinting of fractures, dislocations and sports injuries.
- The product is mainly intended to be used in emergency care in hospitals and during first aid at the scene of accidents.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- For traceability purposes, the product may only be stored in its original packaging. If this is not possible for the user, the user is responsible for correct labelling and traceability.
- There should be sufficient experience in the use of the product in order to avoid injury to patients.
- Do not use the product in case of fractures of the spine or pelvis.
- Do not use the product on open wounds! Danger of infection!
- When cutting the splint to size, sharp edges may occur which must be either covered or bended over!
- For longer periods of application (> 1 hour), a padding underlay should be used (e.g. cloth, dressings or bandages).

3 Modification of the product

- Apart from the described maintenance instructions, the product must not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- The product must be stored away from sunlight, dust, moisture and contamination at temperatures between -25°C and $+65^{\circ}\text{C}$.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ **Damaged products must not be used!**
They must be disposed of in accordance with the instructions in chapter 9.
- The product contains PE and Aluminium.
- The product is radiolucent.
- Please note that the product is supplied in a non-sterile condition. The hygienic preparation of the product lies within the clinical responsibility of the user! The plastic product can be disinfected with a suitable cleaning agent (e.g. alcohol) before use. Sterilization is not permitted!

5 Using the product

- Shape the splint to ensure an optimal fit to the body part.
- Remove any clothing from the affected part of the body.
- Place the splint in the desired location and secure it.
- For longer periods of use, padding on the patient side is recommended.
- For examples of use, see the graphic on the cover page. Other applications are possible.
- When the splint is fixed in place, the affected limb must be exposed in order to observe peripheral blood flow.

6 Product Maintenance/Reprocessing

- The product is reusable subject to risk assessment by the qualified user.
- The service life of reusable products is generally determined by wear and tear and damage caused by use of the product. For this reason, a visual inspection must be carried out before each use, to determine possible material fatigue or impaired functionality.
- The processing of the product lies within the clinical responsibility of the user.
- The product is suitable for manual wipe disinfection.
- For this purpose germicidal cleaning cloths such as Bomix Plus or Korsolex Extra of the X-Wipe-System can be used.
- Then wipe the product with clear water and a clean, soft cloth.
- Sterilization is not permitted!

7 Accessories/Package content

- The user is responsible for the correct combination with other medical products.
- In case of compatibility questions please contact customer service.

8 Replacement parts and repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.
- In the event of a return for repair or replacement the product must be cleaned beforehand in order to exclude any risk to the manufacturer's employees. We reserve the right to reject soiled or contaminated products for safety reasons.

9 Disposal

- Safe disposal is carried out in accordance with current national regulations.

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