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REV GA/2024-06

Besmed Health Business Corp.
No. 5, Lane 116, Wu-Kong 2nd Road,
Wu-Ku District
24888 New Taipei City
Taiwan (R.O.C.)

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GB Contents

1	General	Page	18
2	Safety and Responsibilities	Page	18
3	Modification of the product	Page	19
4	Before using the product	Page	19
5	Using the product	Page	19
6	Product maintenance/Reprocessing	Page	20
6.1	Warnings	Page	20
6.2	Restrictions regarding reprocessing	Page	20
6.3	Preparations at the place of use	Page	21
6.4	Cleaning preparation	Page	21
6.5	Manual cleaning and disinfection	Page	21
6.6	Automatic cleaning and disinfection	Page	22
6.7	Drying	Page	23
6.8	Maintenance, control and inspection	Page	23
6.9	Packaging	Page	23
6.10	Sterilisation	Page	23
6.11	Additional Information	Page	24
7	Accessories/Package content	Page	24
8	Replacement parts and repairs	Page	24
9	Disposal	Page	24
10	Explanation of symbols	Page	25

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1 General

- The product may only be used by trained professionals in cardiopulmonary reanimation and artificial ventilation.
- The anaesthetic masks are intended to be used in combination with a ventilator to ventilate a patient.
- The product is mainly used in hospitals and in emergency care.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- The user should have sufficient experience and knowledge in the assembly, disassembly, cleaning and the use of the product.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- Keep the product free of oil and grease, as this could present an explosion hazard. This includes the exposure to cleaning alcohols, hand cream and adhesives – in contact with highly compressed oxygen explosive reactions may occur.
- Only use replacement parts and accessories provided by Droh. Otherwise the performance of the products may be reduced.
- The reusable products are delivered non-sterile and must therefore be unpacked and fully reprocessed (according to chapter 6) before use.
- The processing instructions under chapter 6 must be followed to avoid any exposure to danger! In case of questions please contact customer service.
- When used in patients with a suspected prion disease, there may be a high risk of infection. It lies within the discretion of the physician to properly dispose of the product according to chapter 9 or to reprocess it according to national regulations (see chapter 6).

3 Modification of the product

- Apart from the described maintenance instructions, the product must not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- Do not expose the product to direct sunlight or other light sources. Store and transport the product in its original packaging at -40°C to $+70^{\circ}\text{C}$.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ **Damaged products must not be used!**
After consulting customer service, the product must be disposed of according to chapter 9.
- The product contains silicone and polysulfone.
- The product is delivered in a non-sterile condition. Cleaning and disinfection must be performed before each use according to instructions chapter 6.
- Attach the mask to a suitable ventilator.
- The product is now ready for use.

5 Using the product

- The operating temperature should be between -40°C and $+50^{\circ}\text{C}$.
- Place the mask over the patient's face, so that no air can escape from the mask.
- Depending on the model a cannula (Luer connector) can be used to adjust the inflation of the cushioned rim.
- Choose the correct type and size of mask according to the anatomical features of the patient.
If applicable observe national legislation.

6 Product maintenance/Reprocessing

- The product is reusable subject to risk assessment by the qualified user.
- The service life of reusable products is generally determined by wear and tear and damage caused by use of the product.
- However, the service life is limited to a maximum of two years or 30 reprocessing cycles from the first usage onwards. The actual wear and tear and the associated number of reprocessing cycles will vary.
For this reason, a functional test must be carried out before each use to determine possible material fatigue or impaired functionality.
- Observe the following procedures for reprocessing. These are validated by the manufacturer.
- A complete reprocessing cycle always consists of cleaning and disinfection (manual or mechanical) and sterilization.

6.1 Warnings

- The product is delivered non-sterile.
- The product must be processed before initial use and before each repeat use.
→ **The delivered packaging cannot be sterilized!**
- The products may only be processed by specifically trained personnel.
- Do not use any pointed or sharp objects as well as hard brushes, scouring agents or other materials that could damage the surface materials!
- The manufacturer accepts no responsibility for damage due to incorrect processing.
- The processing of the product lies within the clinical responsibility of the user.
- No residues may remain in or on the product.
- Do not use damaged products. If necessary, contact customer service.

6.2 Restrictions applicable to reprocessing

- The use of oils and greases is not recommended for silicone products.
- Do not use disinfectants based on formaldehyde, glutaraldehyde, glucoprotamine, cocodipropylenetriamine, quaternary ammonium compounds, phenol derivatives or iodine-cleaving compounds.
- In addition, it must be ensured that no hard or sharp-edged objects come into contact with the individual parts during disinfection and that the parts are rinsed thoroughly afterwards.

6.3 Preparations at the place of use

- It is recommended to rinse the product thoroughly with warm water and a diluted sodium bicarbonate solution (8-10% by volume) immediately after use until all visible foreign matter has been removed.
- Do not let any residues become dried-on.
- Immediately after use, manual disinfection should be performed to reduce the risk of infection for users.

6.4 Cleaning Preparation

- Disconnect the mask from the ventilation system.
- For models with a ventilation plug, remove or open it before machine cleaning, disinfection and sterilisation.
- In the case of models with a ventilation blocker, this must be attached/inserted.
- Remove the cushion and the dome from mask.

6.5 Manual cleaning and disinfection

- Solely manual reprocessing is not recommended due to the lower effectiveness and reproducibility. It should only be used if no automated procedure is available.
- Use fully demineralised water.
- Cleaning and disinfecting solution must be prepared before each manual cleaning cycle. Please follow the manufacturer's instructions for the solution.
- A solution of Sekusept Aktiv is recommended.
- Place the product in the cleaning and disinfecting solution. Pay attention to the manufacturer's recommended concentration and soaking time.
- Even when cleaning manually, make sure that all cavities are properly rinsed through.
- A small brush with soft bristles or a soft sponge can be used to help clean the product.
- After cleaning, rinse sufficiently with fully demineralised water to remove any residues.
- Avoid any water accumulation.
- Dry the product completely.
- Check visually for residues. Repeat the rinsing process if necessary.

6.6 Automatic cleaning and disinfection

- Use fully demineralised water.
- If using alkaline cleaning agents, neutralisation must be carried out first.
- Do not use any desiccants.
- It is recommended to use an automated/ mechanical cleaning and disinfection procedure.
- The product is suitable for machine processing/thermal disinfection ($\leq 93^{\circ}\text{C}$).
- It must be ensured that the selected washer-disinfector complies with ISO 15883.
- Follow the operating instructions of the respective device model.
- Use cleaning agents that are suitable for stainless steel and polysulfone.
- Recommended cleaning agents:
Cleaner: Sekumatic ProClean // Neutraliser: ECOLAB Sekumatic FNZ
- The operating cycle (cleaning, disinfection, rinsing, drying if necessary) must achieve the specified performance.
- Do not overload the basket and avoid shadowing.
- Load products with openings and recesses with open side facing downwards to avoid accumulation of cleaning water.

Instructions:

1. Place the products in the washer-disinfector / arrange around the spiral rack (ensure that all cavities are effectively flushed through).

Observe the loading pattern of your washer-disinfector if necessary.

2. Start a cycle with the following parameters: (e.g. VARIO TD)

- | | |
|----------------------------------|------------------------|
| - Pre-rinse with purified water | 20°C // 1–4 min. |
| - Cleaning with alkaline cleaner | max. 55°C // 5–10 min. |
| - Neutralisation | 20°C // 2–6 min. |
| - Intermediate rinsing | 20°C // 2–3 min |
| - Disinfection | 93°C // 5 min. |
| - Drying | max. 100°C // 30 min |

3. Inspect for any residual soil. Repeat the cycle if necessary.

4. Inspect the products and prepare for sterilization

6.7 Drying

- After cleaning, the products must be dried thoroughly (preferably supported by compressed air).
- Do not use drying agents or cloths which could leave residues on the silicone.

6.8 Maintenance, control and inspection

- Inspect the products after disinfection and before each use for intactness, leakage and cleanliness.
- Check the product for cracks and material changes.
- If residual contamination is still visible, repeat the entire procedure until no more residues are visible.
- Re-attach the cushion to the mask. Pull the nipples completely through the openings.

01-3M

01-4M

01-5M

- When inserting the dome, make sure that the thick edge protrudes from the top of the mask dome.
- After reprocessing and before each use, perform a functional test and check the product for damage.
- For masks with a ventilation plug, close the plug. Now the cushion must not deflate.

6.9 Packaging

- After processing, package the product according to ISO 7153.
- Store sterilised products in sterile packaging in a closed cabinet, protected from dust, light, moisture and temperature fluctuations.

6.10 Sterilisation

- The product must be sterilised before each use.
- Before sterilisation, pack the products in a suitable container or in suitable sterilisation wrapping.
- The correct size of container or wrapping must be chosen.
- Use a sterilisation indicator for the wrapping and note the sterilisation and expiry date on the wrapping.

- The product is only suitable for steam sterilization in a fractionated vacuum process.
Maximum exposure time: 5 min.
Maximum temperature: 134°C
- Follow the instructions of the autoclave manufacturer carefully.
- After steam sterilisation, the product must be allowed to cool down at room temperature.
- The manufacturer accepts no liability for risks arising from improper sterilisation (contamination).

6.11 Additional Information

- Frequent sterilisation can accelerate the natural ageing of the product. In addition, the service life of sterilised products depends very much on carefully performed pre-treatment (cleaning and disinfection).
- When sterilising several products in one sterilisation cycle, the maximum load of the steriliser must not be exceeded.

7 Accessories/Package content

- In case of compatibility questions please contact customer service.


8 Replacement parts and repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.
- If the cause of the error cannot be eliminated, the product must not be put into operation. After consulting the customer service, the product must be disposed of according to chapter 9 or returned to the sales distributor.
- In the event of a return for repair or replacement the entire reprocessing process must be carried out beforehand in order to exclude any risk to the manufacturer's employees.
We reserve the right to reject soiled or contaminated products for safety reasons.

9 Disposal

- If any of the above-mentioned errors have occurred, the product must be disposed of.
- The product must be disposed of at the latest **two years after first use or following 30 reprocessing cycles.**
- Disposal is carried out in accordance with current national regulations.

10 Symbolerklärung / Objaśnienie symboli / Explanation of symbols

		
 Deutschland Kennzeichen	Symbol Niemiec	Symbol for Germany
 Polen Kennzeichen	Symbol Polski	Symbol for Poland
 Groß-Britannien Kennzeichen	Symbol Wielkiej Brytanii	Symbol for Great Britain
 Gebrauchsanweisung beachten	Przestrzegać instrukcji użytkowania	Observe the operating instructions
 Warnhinweis	Wskazówka ostrzegawcza	Warning
 Hersteller	Producent	Manufacturer
 EU-Repräsentant	Przedstawiciel UE	EU Representative
 CE-Symbol	Symbol CE	CE symbol
 Kein Öl oder Schmiermittel verwenden	Nie oliwić	Do not use any grease or other form of lubricant
 Medizinprodukt	Wyrób medyczny	Medical Device
 Importeur	Importer	Importer