Wolfram Droh GmbH, Dammweg 11, 55130 Mainz, Germany

WDI med Sp. z o.o.

ul. Pastelowa 3

70-704 Szczecin

POLSKA


Tel. + 48 91 463 05 71 • Fax + 48 91 463 04 78

info@wdimed.pl • www.wdimed.pl

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(GB) Speed Clip belt sets

REV GA / 2022-07

Anping Guardian Medical Equipment Co., Ltd.
500 meters, Nan Sucun, Anping County
Hengshui City, Hebei Province, 053600
P. R. China

MedNet EC-REP GmbH
Borkstrasse 10
48163 Muenster
Germany

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1 General

- The product may only be used by trained and qualified personnel.
- The product is intended to be used for the complete immobilization of the spine or spinal cord in case of injury.
- The product is intended to be used for secure fixation of a patient on an immobilization or transport device, such as a spineboard.
- The product is mainly used in emergency care and during first aid at the scene of accidents.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- For traceability purposes, the product may only be stored in its original packaging. If this is not possible for the user, the user is responsible for correct labelling and traceability.
- There should be sufficient experience in the use of the product in order to avoid injury to patients.
- Contraindications are not known.
- Worn and damaged belts must be replaced immediately to avoid injury to the patient.
- Make sure that the correct size is used for each individual patient.

3 Modification of the product

- Apart from the described maintenance instructions, the product must not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.


4 Before using the product

- Do not expose the product to direct sunlight or other light sources. They must be stored and transported in a dry place and in their original packaging.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ **Damaged products must not be used!**
They must be disposed of in accordance with the instructions in chapter 9.
- Hooks and fastener are made of chrome-plated steel and plastic. The belt contains nylon and polyester.
- Please note that the product is supplied in a non-sterile condition.
- The hygienic preparation of the product lies within the clinical responsibility of the user! The plastic product can be disinfected with a suitable cleaning agent (e.g. alcohol) before use. Sterilization is not permitted!

5 Using the product

- Immobilize the cervical spine with an appropriate cervical immobilization device.
- With a sufficient number of helpers, first turn the patient onto their side and check for back injuries.
Remove any items from the patient's back pockets.
- Place the spine board against the patient's back.
Roll the patient and spine board back at the same time.
- Position the patient centrally on the board with a sufficient number of helpers. Use the correct method of gripping the patient under the armpits.
- Secure the patient to the spine board with straps or a harness to prevent the patient from moving in any direction (up, down, left or right). Start at the chest area and work your way down. Do not fasten the shoulder straps until last.
- Attach the belts to the pins inside the handle recesses of the spineboard using the speed clips.
- The belts shall be placed over the chest, on the hips, above or beneath the knees of the patient.
- It is advisable to fixate the arms of unconscious patients.
- The positioning may vary caused by individual injuries or local medical guide lines.

- Secure the belts and make sure that the breathing of the patient is not affected. Follow the current medical guide lines.
- Immobilize the head with a head immobilizer on the spine board.
- The head, neck, upper body and pelvis must be immobilized in neutral alignment.

- To close the belt, insert the clip at one end into the clip opening at the other end. If you hear a "click", the belt is closed properly.
- To open the strap, press the button labeled "Press" and release the strap from the buckle.
- To adjust the strap length, adjust the loose end of the belt. To tighten, pull the loose end. To loosen, place the metal buckle at 90° to the taut part of the belt and pull the metal buckle in the opposite direction to the taut end of the belt until the desired length is achieved.
Maximum length: 170 cm
Minimum length: 55 cm

6 Product Maintenance/Reprocessing

- The product is reusable subject to risk assessment by the qualified user.
- The service life of reusable products is generally determined by wear and tear and damage caused by use of the product. For this reason, a visual inspection must be carried out before each use, to determine possible material fatigue or impaired functionality.
- The processing of the product lies within the clinical responsibility of the user.
- The manufacturer accepts no responsibility for damage due to incorrect processing.
- Place the product in the cleaning and disinfecting solution. Pay attention to the manufacturer's recommended concentration and soaking time.
- A solution of Sekusept Plus is recommended.
- A small brush with soft bristles or a soft sponge can be used to help clean the product.
- Check visually for residues. Repeat the rinsing process if necessary.
- Alternatively, clean the product surface with a cloth dampened in a mild soap solution (1:16 ratio) or a solution of household bleach and water (1:100 ratio).
- Do not use aggressive cleaning agents.
- After cleaning, rinse sufficiently with fully demineralised water to remove any residues.
- Sterilization is not permitted!
- Dry the product completely.

7 Accessories/Package content

- The product can be combined with an immobilization and transport device for patient fixation, as well as with a head immobilizer.
- The user is responsible for the correct combination with other medical products.

8 Replacement parts and repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.

9 Disposal

- Safe disposal is carried out in accordance with current national regulations.

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