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
REV 2023-05

Huakun Medical Equipment Co., Ltd.
No. 208, Binjian Road
Cixi High-tech, Ind. Dev. Zone
315301 Zhejiang
CHINA

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Contents

1	General	Page 23
2	Safety and Responsibilities	Page 23
3	Modification of the product	Page 24
4	Before using the product	Page 24
4.1	Correct function of the inlet valve during inspiration	Page 25
4.2	Correct function of the inlet valve during expiration	Page 25
4.3	Correct function of the oxygen connection nozzle	Page 25
4.4	Tightness/correct assembly of the patient valve	Page 26
4.5	Correct function of the pressure relief valve	Page 26
4.6	Correct function of the oxygen reservoir bag (not included)	Page 26
5	Using the product	Page 27
5.1	PVC breathing bag	Page 28
5.2	Ventilation valve with pressure limiting valve	Page 28
5.3	Oxygen supply and oxygen reservoir system	Page 29
5.4	PEEP-valve (not included)	Page 29
6	Product Maintenance / Reprocessing	Page 29
7	Accessories/Package content	Page 30
8	Replacement parts and repairs	Page 30
9	Disposal	Page 30
10	Technical Specifications	Page 31
11	Explanation of symbols	Page 32

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1 General

- The product may only be used by trained professionals in cardiopulmonary reanimation and artificial ventilation.
- The ventilator bags are easy-to-use and functionally reliable hand-held ventilators for the ventilation of patients via a mask, tracheal tube or other ventilation aid (select appropriate model according to patient size).
- Atmospheric air as breathing gas used in ventilation can be augmented with oxygen through the use of the oxygen reservoir system. For this purpose, the system can be connected to an oxygen supply via a standard oxygen tube.
- The product is mainly used in hospitals and in emergency care.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- The user should have sufficient experience and knowledge in the assembly, disassembly and the use of the product.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- The product draws air from the immediate surrounding and must therefore not be used in a toxic atmosphere.
- Keep the product free of oil and grease, as this could present an explosion hazard. This includes the exposure to cleaning alcohols, hand cream and adhesives – in contact with highly compressed oxygen explosive reactions may occur.
- Ensure medically correct use of the pressure limiting valve.
Ventilation pressure too high can result in injuries, including to the lungs.

- Only use the accessories listed under Chapter 7 if expiratory/inspiratory resistance is to be increased.
- Only use replacement parts and accessories provided by Droh. Otherwise, the performance of the products may be limited.
- Always observe the movement of the patient's chest during ventilation. The sound of expiration must always be heard at the patient valve.
- In case of a malfunction of the breathing bag, switch immediately to mouth-to-mouth resuscitation to avoid brain damage to the patient.
- **Single use products may not be reused!**
They must be disposed of immediately after use in accordance with chapter 9.
 The reuse of single use products represents a potential risk for the patient or the user. It can result in contamination and/or reduced functionality.
- Contamination and/or limited functionality of the product can lead to injury, illness or death of the patient.

3 Modification of the product

- Apart from the described maintenance instructions, the product must not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- Do not expose the product to direct sunlight or other light sources. Store and transport the product in its original packaging at -20°C to $+50^{\circ}\text{C}$.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Verify that parts or components do not show any damage, breakage or cracks. The product must also not be contaminated or soiled.
 If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.

→ Damaged products must not be used!

After consulting customer service, the product must be disposed of according to chapter 9.

- Please note that the product is supplied in a non-sterile condition.
- The hygienic preparation of the product lies within the clinical responsibility of the user! The plastic product can be disinfected with a suitable cleaning agent (e.g. alcohol) before use.
Sterilization is not permitted!
- Connect the supplied mask to the patient valve.
- Connect the supplied reservoir bag to the inlet valve.
- When completely assembled, check that all parts are connected tightly.
- The product is now ready for use.
- In addition, perform regular function checks as described below.
- All function checks, unless stated otherwise, shall be performed without the mask and reservoir system.

4.1 Correct function of the inlet valve during inspiration:

- Use the breathing bag with the inlet valve opened.
- Deactivate the pressure limiting valve according to chapter 5.2.
- Squeeze the bag with one hand: The air will escape.
- Now hold the patient valve of the bag closed with your other hand and release your hand holding the bag compressed.
- The bag should now rapidly fill up with air once more.

4.2 Correct function of the inlet valve during expiration:

- Use the breathing bag with the inlet valve opened.
- Deactivate the pressure limiting valve according to chapter 5.2.
- Now hold the patient valve closed with one hand and try to squeeze the bag with your other hand.
- If the bag is not able to be squeezed, then no air can escape via the inlet valve.

4.3 Correct function of the oxygen connection nozzle:

- Deactivate the pressure limiting valve according to chapter 5.2.
- Squeeze the breathing bag together and then seal the reservoir valve with your other hand.
- Release the hand you used to hold the bag compressed.
- If the bag now slowly re-inflates, the oxygen connection is open.

4.4 Tightness/correct assembly of the patient valve:

- Connect a test bag to the patient side of the patient valve (e.g. a phantom lung or anaesthetic mask, choosing the correct size for the size of the breathing bag), and hold the test bag connected to the patient valve in such a way that it is airtight.
- Squeeze the breathing bag several times and observe what the test bag feels like.
- Now press the inflated test bag together and make sure that the air from the test bag escapes via the disc membrane on the patient valve and does not flow back into the breathing bag.

4.5 Correct function of the pressure relief valve:

- Activate the pressure limiting valve according to chapter 5.2.
- Now seal the patient side of the patient valve with your thumb and squeeze the breathing bag several times.
- The pressure relief valve must be heard to function on every expiration.

4.6 Correct function of the oxygen reservoir bag:

- Connect the reservoir system via the reservoir valve to the breathing bag (see also the description under chapter 5.3).
- Connect the breathing bag including the reservoir system to an oxygen supply via the oxygen valve.
- Open the oxygen supply to a rate of 15 l/min.
- With one hand seal the outlet of the patient valve.
- Make sure that the reservoir bag inflates and deflates again when the breathing bag is compressed.

5 Using the product

- The selected model must be suitable for the patient.
For the available models see chapter 10.
- Select the ventilation technique, ventilation frequency and, if necessary, ventilation pressure according to the current applicable guidelines.
- The following is an example procedure for ventilating an adult patient:

– First remove any foreign objects from the patient's mouth and throat.

– Turn the patient onto his/her back and position yourself behind the head. Hold the patient's head stretched backwards by pulling the chin towards you. This will open the patient's airways.

– Press the mask firmly over the patient's face using your thumb and index finger. At the same time use your other fingers to support the chin.
– Squeeze the bag with your other hand to ventilate the patient.

- After the compression phase, the elastic breathing bag unfolds automatically and inflates once more via the inlet membrane of the inlet valve that opens as a result of the vacuum created.
- Please observe:
 - a) whether the patient's chest rises and falls in step with your ventilation actions.
 - b) the colour of the patient's lips and face through the transparent face mask.
 - c) whether the mask mists over during exhalation.
 - d) whether the patient valve works correctly.
- In case the patient vomits during ventilation, remove the vomit and clean the patient's airways immediately. Before continuing with the ventilation, squeeze the bag a few times.
- A patient with endotracheal intubation can be connected directly to the outlet on the patient valve.

5.1 PVC breathing bag

- A quantity of gas (air or oxygen-enriched air mixture) flows through the ventilation valve to the patient when the bag is manually compressed.
- After the compression phase, the elastic breathing bag unfolds automatically and inflates once more via the inlet membrane of the inlet valve that opens as a result of the vacuum created.
- Expiratory membranes prevent expiratory air from flowing back into the breathing bag.
- When the patient breathes independently, the inlet valve also opens as a result of the vacuum caused in the bag, allowing the patient to inhale.

5.2 Ventilation valve with pressure-limiting valve


- The pressure-limiting valve allows ventilation pressure to be reduced to a maximum of 60 mbar for the adult model and 40 mbar for the child and infant model.

The valve will open automatically when these pressures are exceeded.

- The valve is set to activated by default.
- To deactivate, carefully press the pressure limiting valve inwards. Turn the valve by 90° (the valve can be turned in a clockwise or anti-clockwise direction).
- The correct setting of the pressure limiting valve can be checked as described in chapter 4.5.

5.3 Oxygen supply and oxygen reservoir system

- When additional oxygen is administered, the oxygen supply tube is connected to the small nozzle on the inlet valve.
- A reservoir system is required in order to achieve high oxygen content and to avoid any oxygen loss.
- The oxygen reservoir system consists of a soft bag (reservoir bag) that is only connected to the reservoir valve.
The complete system is then connected via the adapter to the large cone of the inlet valve.
- The specifications for the oxygen concentration are listed in the table in chapter 10.

5.4 PEEP-valve (not included)

- The breathing bag can be combined with a PEEP valve.
In the event of use follow the corresponding instructions for use.

6 Product maintenance/Reprocessing

- Single use products may not be reused!
- After termination of the application, the product must be disposed of according to chapter 9.

7 Accessories/Package content

- The following three models are available. The selection is according to the weight of the patient. Observe the overview given in Chapter 10.
- All models comprise a mask, a patient valve, self-inflating PVC breathing bag, inlet valve and reservoir bag.
- A PEEP valve can be connected to all models in combination with an adapter.
- Only use replacement parts and accessories provided by Droh.
- In case of compatibility questions please contact customer service.

8 Replacement parts and repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.
- If the cause of the error cannot be eliminated, the product must not be put into operation. After consulting the customer service, the product must be disposed of according to chapter 9 or returned to the sales distributor.

9 Disposal

- The product must be disposed of after use.
- Please note the expiry date if applicable.
The product must not be used after the expiry date!
- Safe disposal is carried out in accordance with current national regulations

10 Technical Specifications

- The data recorded and listed as part of the specifications are mean values and may vary owing to changed conditions (e.g. temperature fluctuations/different levels of strength exerted by users when squeezing the breathing bag etc.).

	Infant breathing bags	Child breathing bags	Adult breathing bags
Patient weight	up to 10 kg	8 kg to 23 kg	over 23 kg
Bag Weight	350 g	410 g	600 g
Volume	310 ± 80 ml	650 ± 100 ml	1600 ± 150 ml
Dimensions (L x W x H)	220 x 100 x 100 mm	260 x 110 x 110 mm	330 x 140 x 140 mm
Reservoir volume	1600 ± 200 ml	1600 ± 200 ml	2000 ± 200 ml
Reservoir volume	30 - 45 cm H ₂ O (15 L/min)	30 - 45 cm H ₂ O (15 L/min) 40 - 60 cm H ₂ O (60 L/min)	40 - 60 cm H ₂ O (60 L/min)
Max. tidal volume /ml	150 ml	350 ml	600 ml
Inspiratory resistance	≤ 5 cm H ₂ O (5 L/min)	≤ 5 cm H ₂ O (50 L/min)	≤ 5 cm H ₂ O (50 L/min)
Expiratory resistance	≤ 5 cm H ₂ O (5 L/min)	≤ 5 cm H ₂ O (50 L/min)	≤ 5 cm H ₂ O (50 L/min)
Patient valve dead space	≤ 20 ml	≤ 40 ml	≤ 65 ml
Expiratory port	15 mm conical socket and 22 mm OD conical connector acc. to ISO 5356		
Connection for oxygen supply	6 mm nozzle		
Operating temperature	-18°C to +50°C		

Materials

Rigid plastic parts	Polycarbonate
Pressure limiting valve spring	Stainless steel
Reservoir bag, breathing bag	Ethylene vinyl acetate (EVA)

11 Symbolerklärung / Objaśnienie symboli / Explanation of symbols

		
 Deutschland Kennzeichen	Symbol Niemiec	Symbol for Germany
 Polen Kennzeichen	Symbol Polski	Symbol for Poland
 Groß-Britannien Kennzeichen	Symbol Wielkiej Brytanii	Symbol for Great Britain
 Gebrauchsanweisung beachten	Przestrzegać instrukcji użytkowania	Observe the operating instructions
 Warnhinweis	Wskazówka ostrzegawcza	Warning
 Hersteller	Producent	Manufacturer
 EU-Repräsentant	Przedstawiciel UE	EU Representative
 CE-Symbol	Symbol CE	CE symbol
 Zur einmaligen Verwendung/ nur zur Verwendung an einem Patienten	Do jednorazowego użytku/tylko do użytku na jednym pacjencie	Single-use only/Only to be used on one patient
 Kein Öl oder Schmiermittel verwenden	Nie oliwić	Do not use any grease or other form of lubricant
 Nicht verwenden wenn Packung beschädigt ist	Nie użytkować, jeśli opakowanie jest uszkodzone	Do not use if package is damaged
 Medizinprodukt	Wyrób medyczny	Medical Device
 Importeur	Importer	Importer