Wolfram Droh GmbH, Dammweg 11, 55130 Mainz, Germany

WDI med Sp. z o.o.

ul. Pastelowa 3

70-704 Szczecin

POLSKA

Tel. +48 91 463 05 71 • Fax +48 91 463 04 78

info@wdimed.pl • www.wdimed.pl

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1 General

- The product may only be used by trained and qualified personnel.
- At least two people are required for safe and correct use.
- The product is intended to be used for immobilization of cranial and cervical spine stability during transportation, especially for patients with suspected cranial or spinal trauma.
- The product is usually combined with a scoop stretcher, a spine board or other combi-carriers to ensure the most effective immobilization during transport or rescue.
- The product is mainly intended to be used in emergency care and during first aid at the scene of accidents.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- For traceability purposes, the product may only be stored in its original packaging. If this is not possible for the user, the user is responsible for correct labelling and traceability.
- There should be sufficient experience in the use of the product in order to avoid injury to patients.
- Fire Hazard: The product must not be exposed to heat sources or flammable substances.
- Do not apply excessive force at any time during use to avoid injury.
- The product must not come into contact with pointed or sharp objects to avoid cracks or other damage.
- A fully functioning spare of the product should always be kept ready to hand where the product is being used.
- In case of obvious malfunction of the product, it must not be used! It must be replaced immediately.
- Worn and damaged belts must be replaced immediately to avoid injury to the patient.
- Make sure that the correct size is used for each individual patient.

3 Modification of the product

- Apart from the described maintenance instructions, the product must not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product


- Do not expose the product to direct sunlight or other light sources. They must be stored and transported in a dry place and in their original packaging.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ **Damaged products must not be used!**
They must be disposed of in accordance with the instructions in chapter 9.
- The product contains PE and Vinyl.
- The product is suitable for X-ray and MRI examinations.
- Please note that the product is supplied in a non-sterile condition.
- The hygienic preparation of the product lies within the clinical responsibility of the user! The plastic product can be disinfected with a suitable cleaning agent (e.g. alcohol) before use.
Sterilization is not permitted!
- Simultaneous use of a cervical collar is recommended.

5 Using the product

- Attach the base plate to the spine board.
Top strap to central opening (optional, if available).
Lateral fixation (using included positioning straps).

- Attach the head blocks.
The ears should be positioned within the ear openings.

- Fix the head at the forehead and chin (using included fixation straps).

6 Product Maintenance/Reprocessing

- The product is reusable subject to risk assessment by the qualified user.
- The service life of reusable products is generally determined by wear and tear and damage caused by use of the product. For this reason, a visual inspection must be carried out before each use, to determine possible material fatigue or impaired functionality.
- The processing of the product lies within the clinical responsibility of the user.
- The product is suitable for manual wipe disinfection.
- For this purpose germicidal cleaning cloths such as Bomix Plus or Korsolex Extra of the X-Wipe-System can be used.
- Alternatively, clean the product surface with a cloth dampened in a mild soap solution (1:16 ratio) or a solution of household bleach and water (1:100 ratio).
- Do not use aggressive cleaning agents.
- Then wipe the product with clear water and a clean, soft cloth.
- Sterilization is not permitted!

7 Accessories/Package content

- Base plate (1 x)
- Head block/pad (2 x)
- Positioning strap (2 x)
- Patient fixation strap for forehead/chin (2 x)

8 Replacement parts and repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.
- In the event of a return for repair or replacement the product must be cleaned beforehand in order to exclude any risk to the manufacturer's employees. We reserve the right to reject soiled or contaminated products for safety reasons.

9 Disposal

- Safe disposal is carried out in accordance with current national regulations.

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