Wolfram Droh GmbH

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
MD

REV 2024-09

Hangzhou Formed Medical Devices Co., Ltd.
Building 3, No. 292 Gonghe North Road
311404 Hangzhou
P. R. China

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(GB) Contents

1	General	Page	17
2	Safety and Responsibilities	Page	17
3	Modification of the product	Page	19
4	Before using the product	Page	19
5	Using the product	Page	20
6	Product Maintenance / Reprocessing	Page	21
7	Accessories	Page	21
8	Replacement parts and repairs	Page	22
9	Disposal	Page	22
10	Technical Specifications	Page	22
11	Error handling	Page	23
12	Explanation of symbols	Page	24

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1 General

- The medical products may only be used by trained and qualified personnel.
- Use in neonates and infants must only be performed by anaesthesiologists with sufficient experience.
- The laryngeal mask is an alternative to the simple face mask and is used for routine and critical anaesthesia procedures to secure the airway.
- The product can also be inserted for resuscitation of deeply unconscious patients with absence of glossopharyngeal and laryngeal reflexes to secure the airway.
- Furthermore, it can be used when tracheal intubation must be excluded or has failed.
- The product is mainly used in emergency care for intubation of accident victims or for intubation before surgery in the field of anaesthesia.
- Laryngeal masks consist as a complete system of three parts: tube (with 15 mm standard connector), cuff and tube for filling the cuff. All components are latex-free.
- The mask is shaped to conform to the contours of the laryngopharynx, with the lumen opening into the larynx.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! PLEASE READ CAREFULLY BEFORE USE
PLEASE KEEP FOR FURTHER REFERENCE

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- There should be sufficient experience in the use of the product in order to avoid injury to patients.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.

- The product draws air from the immediate surrounding and must therefore not be used in a toxic atmosphere.
- Use of the product in the presence of oxygen or nitrous oxide, during procedures involving lasers or electrosurgical probes, may cause the product to ignite. This can result in burns and poisoning of the patient.
- The product is not CT/MRI compatible.
- Do not immerse or soak the product in liquid.
- The product is contraindicated in the following situations:
 - unconfirmed fasting
 - morbid obesity
 - massive and/or abdominal or thoracic injuries
 - pregnancy after 14 weeks of pregnancy
 - conditions with delayed gastric emptying or opioid medication prior to fasting
 - persistent reduced lung compliance (the laryngeal mask forms a low-pressure seal around the larynx)
- The cuff must be completely deflated during intubation and extubation. Otherwise, it can lead to serious injuries to the respiratory tract.
- If the laryngeal mask is used in deeply unconscious patients requiring resuscitation, there is a risk of regurgitation or aspiration, which must be weighed against the benefit of securing the airway.
- If ventilation is inadequate, neonates and children will develop oxygen deficiency more quickly due to increased oxygen consumption.
- Laryngeal mask may cause minor (e.g., sore throat) to severe (e.g., aspiration) side effects.
- Take care that the tip of the cuff does not enter the vallecula or the glottis and does not press on the epiglottis or arytenoid cartilage.
- The cuff should be moved against the patient's pharynx during insertion. Resistance can be felt when the mask is fully inserted.
- Do not apply excessive force at any time during intubation to avoid injury.
- A fully functioning spare of the product should always be kept ready to hand where the product is being used.
- Do not use silicone-based lubricants.
- Remove all of the patient's removable dentures before use.
- Once placed in the patient, the mask must not be over-inflated, otherwise it may cause internal injuries and suffocation.

- The product must not come into contact with pointed or sharp objects to avoid cracks or other damage.
- If signs of leakage from the product occur during ventilation, replace the mask immediately.
- Long-term ventilation may lead to damage to the mucous membranes.
- Single use products may not be reused!
- They must be disposed of immediately after use in accordance with chapter 9.
- The reuse of single use products represents a potential risk for the patient or the user. It can result in contamination and/or reduced functionality.
- Contamination and/or limited functionality of the product can lead to injury, illness or death of the patient.

3 Modification of the product

- Apart from the described maintenance instructions, the product may not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- There are no special conditions to be observed during transport.
- The product must be stored away from sunlight, dust, moisture and contamination at temperatures between 10°C and 25°C.
- Make sure that the correct size is used for each individual patient.
- Always observe the maximum cuff volume. Various sizes are available for use in paediatrics.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- A protective double shell made of transparent plastic may be found around the cuff, which serves only as supportive packaging material for transportation. This must not be inserted into the patient's mouth and must be disposed of before use.
- Check the tight fit of the ISO connector.

- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ **Damaged products must not be used!**
They must be disposed of in accordance with the instructions in chapter 9.
- The product contains silicone, PVC, PC, PP and/or ABS.
- Wire-reinforced laryngeal masks also contain stainless steel 304.
- The product is supplied in a sterile condition.
- To minimize contamination, gloves should be worn before and during intubation.
- A function check must be carried out before each use.
If this check fails, the product must not be used.
- Check that the cuff and the inside of the tube are free of obstructions.
Remove any foreign objects from the product.
- Test the cuff and valve by a seal test before use.
 - Use a syringe to vent the cuff via the valve so that its walls are close together.
 - Fill the cuff with approximately 50% more air than the recommended maximum volume.
 - The resulting balloon should be elliptical, not round. There must be no uneven bulges.
 - If there are signs of leakage, do not use the product!
Contact customer service.
 - Then empty the cuff again so that it is completely free of wrinkles.
- The product is now ready for use.

5 Using the product

- During intubation and extubation as well as during ventilation, observe the currently valid medical guidelines and the information in the instructions for use.
- Adequate depth of anaesthesia of the patient must be ensured.
- The product is suitable for spontaneously breathing patients when used with volatile anaesthetics or intravenous narcotics. (Provided the depth of anaesthesia is appropriate for the respective surgical stimulus and the cuff is not over-inflated).
- Note the position of the patient's head. The so-called "sniffing position" is common, as used for tracheal intubation.

- Position the laryngeal mask according to recognized medical techniques such as the index finger insertion technique.
- Insert the product orally into the patient's throat according to the valid procedures. The black X-ray line of the tube should point in the direction of the patient's nose.
- If necessary you can use a water-soluble lubricant. Apply this sufficiently to the back of the tip of the cuff. Follow the manufacturer's quantity instructions and avoid blockage of the lumen and airway.
- After intubation, inflate the cuff with just enough air to obtain a seal.
- The maximum recommended cuff pressure is 60 cmH₂O.
Cuff pressure should be monitored continuously ideally with a cuff pressure gauge. If a cuff-pressure gauge is not available, first fill the cuff with only two-thirds of the recommended maximum volume and then in small increments until you can no longer hear air escaping in the mouth and throat.
- Do not hold the cuff during inflation, as this will prevent the mask from sliding into the correct position.
- During the first three to four breaths, slight leakage may occur.
The mask should then have positioned itself correctly.
If the leak persists, the first step should be to check for low ventilation pressure and adequate depth of anaesthesia.
If in doubt, replace the product.
- The correct position of the mask can be observed as follows:
 - Slight upward movement of the tube when inflating the cuff.
 - A slight oval swelling in the neck in the region of the cricoid and thyroid cartilage, the cuff must not be visible in the oral cavity.
- If repositioning of the mask is necessary, deflate the cuff completely and inflate again after repositioning.
- Do not completely deflate the cuff before extubation to prevent secretions from entering the laryngeal mask.
Perform extubation according to accepted medical guidelines.

6 Product maintenance / Reprocessing

- Single use products may not be reused!
- After termination of the application, the product must be disposed of according to chapter 9.

7 Accessories

- The user is responsible for the correct combination with other medical products.
- In case of compatibility questions please contact customer service.

8 Replacement parts and repairs

- In case of functional or material problems, contact customer service immediately.

9 Disposal

- The product must be disposed of after use.
- Please note the expiry date if applicable.
The product must not be used after the expiry date!
- Safe disposal is carried out in accordance with current national regulations.

10 Technical Specifications

- When selecting the correct mask size, please note that the following specifications are only approximate. The anatomy of the patient shall also be taken into account.
- The volume specifications are maximum test volumes.
These values must not be reached during use on patients.
The inflation volume depends on the seal pressure (see Chapter. 5).

Mask size	1	1,5	2	2,5	3	4	5
Patient weight † (kg)	< 5	5 - 10	10 - 20	20 - 30	30 - 50	50 - 70	> 70
Tube inner diameter (mm)	5,0	6,0	7,0	8,5	10,0	10,0	11,0
Tube outer diameter (mm)	8,7	10,0	11,3	13,3	15,3	15,3	16,7
Max. cuff volume (ml)	4	7	10	14	20	30	40
Cuff inflation valve	Luer connection (DIN 13090 F6)						
Airway connector	15 mm male ISO						

11 Error handling

Problem	Cause	Solution
Leakage beyond the first 3–4 breaths.	Incorrect positioning of the mask.	Remove and reinsert the mask.
Unexpected regurgitation; oxygen saturation within acceptable range.	The depth of anaesthesia is no longer sufficiently guaranteed. An indication of this is the appearance of liquid in the tube. In spontaneously breathing patients, coughing or breath holding may be a sign.	Immediately place patient in head-down position and briefly disconnect mask from anaesthesia machine. Check the depth of anaesthesia and increase intravenously if necessary. Suction the tube.
Signs of leakage during use with IPPV (intermittent positive pressure ventilation).	– Too shallow anaesthesia causing glottis closure. – Inadequate neuromuscular block. – Reduction in lung compliance associated with the surgical/ diagnostic procedure. – Displacement of the mask by turning or pulling.	– Deepen the anaesthesia. – Insufficient relaxation. – Adjustment of the ventilation pressure/PEEP. – Reposition the mask.
If the problem cannot be corrected, the laryngeal mask should be removed and breathing ensured by other means.		

12 Symbolerklärung / objaśnienie symboli / Explanation of symbols

		
 Deutschland Kennzeichen	Symbol Niemiec	Symbol for Germany
 Polen Kennzeichen	Symbol Polski	Symbol for Poland
 Groß-Britannien Kennzeichen	Symbol Wielkiej Brytanii	Symbol for Great Britain
 Gebrauchsanweisung beachten	Przestrzegać instrukcji użytkowania	Observe the operating instructions
 Warnhinweis	Wskazówka ostrzegawcza	Warning
 Hersteller	Producent	Manufacturer
 EU-Repräsentant	Przedstawiciel UE	EU Representative
 CE-Symbol	Symbol CE	CE symbol
 Nicht verwenden wenn Packung beschädigt ist	Nie użytkować, jeśli opakowanie jest uszkodzone	Do not use if package is damaged
 Zur einmaligen Verwendung / nur zur Verwendung an einem Patienten	Do jednorazowego użytku / tylko do użytku na jednym pacjencie	Single-use only / Only to be used on one patient
 Steril durch EO-Begasung	Produkt sterylizowany tlenkiem etylenu	Sterilized by exposure to EO gas
 Nicht erneut sterilisieren	Nie sterylizować ponownie	Do not re-sterilize
 Medizinprodukt	Wyrób medyczny	Medical Device
 Importeur	Importer	Importer