Wolfram Droh GmbH

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1 General

- This product may only be used by trained professionals under medical instruction.
- The product is intended to be used during operations as a drainage aid for the transportation (for max. 29 days) of blood, body fluids or other biocompatible fluids into and out of any human body orifice.
It serves as a connection between a suction catheter and a suction device.
- With Yankauer suction tips, the vacuum can be controlled via the handpiece.
- For suction tubes with fingertip, the vacuum can be controlled by opening and closing the plug.
- It is mainly used in hospitals in the operating room for surgical drainage.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- The product contains di(2-ethylhexyl) phthalate (DEHP) and should therefore only be used temporarily in pregnant women, nursing mothers, babies and young children.
- The product must not come into contact with a pointed or a sharp object which could damage it. This could lead to breakages or similar damage and the tubing would start to leak.
- **Single use products may not be reused!**
They must be disposed of immediately after use in accordance with chapter 9. The reuse of single use products represents a potential risk for the patient or the user. It can result in contamination and/or reduced functionality.
- Contamination and/or limited functionality of the product can lead to injury, illness or death of the patient.

3 Modification of the product

- Apart from the described maintenance instructions, the product may not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.

- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- Do not expose the product to direct sunlight or other light sources. They must be stored and transported in a dry place and in their original packaging.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ **Damaged products must not be used!**
They must be disposed of in accordance with the instructions in chapter 9.
- The product contains PVC.
- The product is supplied sterile.
- After assembly, the product must be checked for sufficient pull-off force and for kinks as well as any cracks or hardening.
- The product is now ready for use.

5 Using the product

- After connection to a suitable terminal device, the product is ready for use.
- Suction tube with vacuum control:
Close the opening at the fingertip outlet with the plug provided for suction. When the opening is not sealed, the tube can be used to control the vacuum.

6 Product maintenance/Reprocessing

- Single use products may not be reused!
- The product must not be used on more than one patient.
- After termination of the application, the product must be disposed of according to chapter 9.

7 Accessories

- The user is responsible for the correct combination with other medical products.

8 Replacement parts and repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.
- If the cause of the error cannot be eliminated, the product must not be put into operation. After consulting the customer service, the product must be disposed of according to chapter 9 or returned to the sales distributor.

9 Disposal

- The product must be disposed of after use.
- Please note the expiry date if applicable.
The product must not be used after the expiry date!
- Safe disposal is carried out in accordance with current national regulations.

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Single-use only
/Only to be used
on one patient

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