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

REV GA 2022-02

HEADSTAR Medical Products Co., Ltd.
9F, No. 288, Sec. 2, New Taipei Blvd.
Xinzhuan Dist., 242 New Taipei City
Taiwan, R.O.C.

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1	General	Page	37
2	Safety and Responsibilities	Page	37
3	Modification of the product	Page	38
4	Before using the product	Page	38
4.1	Correct function of the inlet valve during inspiration	Page	39
4.2	Correct function of the inlet valve during expiration	Page	39
4.3	Correct function of the oxygen connection nozzle	Page	39
4.4	Tightness/correct assembly of the patient valve	Page	39
4.5	Correct function of the pressure relief valve	Page	39
4.6	Correct function of the oxygen reservoir bag (not included)	Page	40
5	Using the product and functional description	Page	40
5.1	Silicone breathing bag	Page	41
5.2	Ventilation valve with pressure limiting valve	Page	41
5.3	Oxygen supply and oxygen reservoir system	Page	42
5.4	PEEP-valve (not included)	Page	42
6	Product maintenance/Reprocessing	Page	43
6.1	Warnings	Page	43
6.2	Restrictions applicable to reprocessing	Page	43
6.3	Preparation at the place of use	Page	44
6.4	Preparation for cleaning	Page	44
6.5	Manual cleaning and disinfection	Page	44
6.6	Automatic cleaning and disinfection	Page	44
6.7	Drying	Page	45
6.8	Maintenance, control and inspection	Page	45
6.9	Packaging	Page	46
6.10	Sterilisation	Page	46
6.11	Additional Information	Page	46
7	Accessories and replacement parts	Page	46/47
8	Repairs	Page	48
9	Disposal	Page	48
10	Technical Specifications	Page	49/50
11	Detailed illustrations	Page	51
12	Explanation of symbols	Page	52

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1 General

- The product may only be used by trained professionals in cardiopulmonary reanimation and artificial ventilation.
- The ventilator bags are easy-to-use and functionally reliable hand-held ventilators for the ventilation of patients via a mask, tracheal tube or ventilation aid (select appropriate model according to patient size).
- Atmospheric air as breathing gas used in ventilation can be augmented with oxygen through the use of the oxygen reservoir system. For this purpose, the system can be connected to an oxygen supply via a standard oxygen tube.
- The product is mainly used in hospitals and in emergency care.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- The user should have sufficient experience and knowledge in the assembly, disassembly, cleaning and the use of the product.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- The product draws air from the immediate surrounding and must therefore not be used in a toxic atmosphere.
- Keep the product free of oil and grease, as this could present an explosion hazard. This includes the exposure to cleaning alcohols, hand cream and adhesives – in contact with highly compressed oxygen explosive reactions may occur.
- Ensure medically correct use of the pressure limiting valve. Ventilation pressure too high can result in injuries, including to the lungs.
- Only use the accessories listed under chapter 7 if expiratory/inspiratory resistance is to be increased.
- Only use replacement parts and accessories provided by Droh. Otherwise the performance of the products may be limited.
- Always observe the movement of the patient's chest during ventilation.

- In case of a malfunction of the breathing bag switch immediately to mouth-to-mouth resuscitation to avoid brain damage to the patient.
- The reusable products are delivered non-sterile and must therefore be unpacked and fully reprocessed (according to chapter 6) before use.
- The processing instructions under chapter 6 must be followed to avoid any exposure to danger! In case of questions please contact customer service.
- Repeated use can lead to hardening and cracking. In this case the product must be disposed of according to chapter 9.
- When used in patients with a suspected prion disease, there may be a high risk of infection.
It lies within the discretion of the physician to properly dispose of the product according to chapter 9 or to reprocess it according to national regulations (see chapter 6).

3 Modification of the product

- Apart from the described maintenance instructions, the product must not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- Do not expose the product to direct sunlight or other light sources.
Store and transport the product in its original packaging at -20°C to $+50^{\circ}\text{C}$.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ Damaged products must not be used!
After consulting customer service, the product must be disposed of according to chapter 9.
- The product contains silicone and polysulfone.
- Connect the supplied mask to the patient valve.
- In case of further accessories (see chapters 7) assemble these according to the detailed illustrations under chapter 11.
- When completely assembled, check that all parts are connected tightly.

- The product is now ready for use.
- In addition, perform regular function checks as described below.
- All function checks, unless stated otherwise, shall be performed without the mask and reservoir system.

4.1 Correct functionality of the inlet valve during inspiration:

- Use the breathing bag with the inlet valve opened.
- Deactivate the pressure limiting valve according to chapter 5.2.
- Squeeze the bag with one hand: The air will escape.
- Now hold the patient valve of the bag closed with your other hand and release your hand holding the bag compressed.
- The bag should now rapidly fill up with air once more.

4.2 Correct function of the inlet valve during expiration:

- Use the breathing bag with the inlet valve opened.
- Deactivate the pressure limiting valve according to chapter 5.2.
- Now hold the patient valve closed with one hand and try to squeeze the bag with your other hand.
- If the bag is not able to be squeezed, then no air can escape via the inlet valve.

4.3 Correct function of the oxygen connection nozzle:

- Deactivate the pressure limiting valve according to chapter 5.2.
- Squeeze the breathing bag together and then seal the reservoir valve with your other hand.
- Release the hand you used to hold the bag compressed.
- If the bag now slowly re-inflates, the oxygen connector is open.

4.4 Tightness/correct assembly of the patient valve:

- Connect a test bag to the patient side of the patient valve (e.g. a phantom lung or anaesthetic mask, choosing the correct size for the size of the breathing bag), and hold the test bag connected to the patient valve in such a way that it is airtight.
- Squeeze the breathing bag several times and observe what the test bag feels like.
- Now press the inflated test bag together and make sure that the air from the test bag escapes via the disc membrane on the patient valve and does not flow back into the breathing bag.

4.5 Correct function of the pressure relief valve:

- This does not apply to models **without** a pressure limiting valve.
- Activate the pressure limiting valve according to chapter 5.2.
- Now seal the patient side of the patient valve with your thumb and squeeze the breathing bag several times.
- The pressure relief valve must be heard to function on every expiration.

4.6 Correct function of the oxygen reservoir bag (not included):

- Connect the reservoir system via the reservoir valve to the breathing bag (see also the description under chapter 5.3).
- Connect the breathing bag including the reservoir system to an oxygen supply via the oxygen valve.
- Open the oxygen supply to a rate of 15 l/min.
- Make sure that the reservoir bag inflates and deflates again when the breathing bag is compressed.

5 Using the product

- The selected model must be suitable for the patient.
For the available models see chapter 10.
- Select the ventilation technique, ventilation frequency and, if necessary, ventilation pressure according to the current applicable guidelines.
- The following is an example procedure for ventilating an adult patient:

– First remove any foreign objects from the patient's mouth and throat.

– Turn the patient onto his/her back and position yourself behind the head. Hold the patient's head stretched backwards by pulling the chin towards you. This will open the patient's airways.

– Press the mask firmly over the patient's face using your thumb and index finger. At the same time use your other fingers to support the chin.
– Squeeze the bag with your other hand to ventilate the patient.

- After the compression phase, the elastic breathing bag unfolds automatically and inflates once more via the inlet membrane of the inlet valve that opens as a result of the vacuum created.
- Please observe:
 - a) whether the patient's chest rises and falls in step with your ventilation actions
 - b) the colour of the patient's lips and face through the transparent face mask
 - c) whether the mask mists over during exhalation
 - d) whether the patient valve works correctly.

- In case the patient vomits during ventilation, remove the vomit and clean the patient's airways immediately. Before continuing with the ventilation, squeeze the bag a few times.
- A patient with endotracheal intubation can be connected direct to the outlet on the patient valve.

5.1 Silicone breathing bag

- A quantity of gas (air or oxygen-enriched air mixture) flows through the ventilation valve to the patient when the bag is manually compressed.
- After the compression phase, the elastic breathing bag unfolds automatically and inflates once more via the inlet membrane of the inlet valve that opens as a result of the vacuum created.
- Expiratory membranes prevent expiratory air from flowing back into the breathing bag.
- When the patient breathes independently, the inlet valve also opens as a result of the vacuum caused in the bag, allowing the patient to inhale.

5.2 Ventilation valve with pressure-limiting valve

- The pressure-limiting valve allows ventilation pressure to be reduced to a maximum of 60 mbar for the adult model and 40 mbar for the child and infant model.
The valve will open automatically when these pressures are exceeded.
- The valve is set to activated by default (see illustration).
- To deactivate, carefully press the pressure limiting valve inwards. Turn the valve by 90° (the valve can be turned in a clockwise or anti-clockwise direction).
- The correct setting of the pressure limiting valve can be checked as described in chapter 4.5.

Pressure-limiting valve closed


Pressure-limiting valve open

5.3 Oxygen supply and oxygen reservoir system

- When additional oxygen is administered, the oxygen supply tube is connected to the small nozzle on the inlet valve.
- A reservoir system is required in order to achieve high oxygen content and to avoid any oxygen loss.
- The oxygen reservoir system consists of a soft bag (reservoir bag) that is only connected to the reservoir valve.
The complete system is then connected via the adapter to the large cone of the inlet valve.
- The specifications for the oxygen concentration are listed in the table in chapter 10.

5.4 PEEP-valve (not included)

- The breathing bag can be combined with a PEEP valve.
In the event of use follow the corresponding instructions for use.

6 Product maintenance / Reprocessing

- The product is reusable subject to risk assessment by the qualified user.
- The service life of reusable products is generally determined by wear and tear and damage caused by use of the product.
- However, the service life is limited to a maximum of **two years or 30 reprocessing cycles** from the first usage onwards.
The actual wear and tear and the associated number of reprocessing cycles will vary. For this reason, a functional test must be carried out before each use to determine possible material fatigue or impaired functionality.
- Observe the following procedures for reprocessing.
These are validated by the manufacturer.
- A complete reprocessing cycle always consists of cleaning and disinfection (manual or mechanical) and sterilization.

6.1 Warnings

- **The reservoir bag is intended for single use only.**
This product must not be reused! Reprocessing is not permitted! It must be disposed of after use according to chapter 9.
- The product is delivered non-sterile.
- The product must be processed before initial use and before each repeat use.
→ **The delivered packaging cannot be sterilized!**
- The products may only be processed by specifically trained personnel.
- Do not use any pointed or sharp objects as well as hard brushes, scouring agents or other materials that could damage the surface materials!
- The manufacturer accepts no responsibility for damage due to incorrect processing.
- The processing of the product lies within the clinical responsibility of the user.
- No residues may remain in or on the product.
- Do not use damaged products.
If necessary, contact customer service.

6.2 Restrictions applicable to reprocessing

- The use of oils and greases is not recommended for silicone products.
- Do not use disinfectants based on formaldehyde, glutaraldehyde, glucoprotamine, cocosdipropylenetriamine, quaternary ammonium compounds, phenol derivatives or iodine-cleaving compounds.
- In addition, it must be ensured that no hard or sharp-edged objects come into contact with the individual parts during disinfection and that the parts are rinsed thoroughly afterwards.
- All parts of the breathing bag are suitable for reprocessing except the reservoir bag (01-19H & 01-21H). This is a PVC product and should be disposed of after use according to chapter 9.

6.3 Preparation at the place of use

- It is recommended to rinse the product thoroughly with warm water and a diluted sodium bicarbonate solution (8-10% by volume) immediately after use until all visible foreign matter has been removed.
- Do not let any residues become dried-on.
- Immediately after use, manual disinfection should be performed to reduce the risk of infection for users.

6.4 Preparation for cleaning

- Disassemble the breathing bag into its individual parts (see detailed illustration under chapter 12).
- Separate the silicone cushion from the dome of the mask.
Do not remove the fastening nipple from the mask!
- Disassemble the inlet valve, patient valve (breathing valve) and reservoir valve.
- Unscrew the pressure limiting valve from the patient valve.
Do not disassemble the pressure limiting valve (01-29H/ 01-31H) itself!
If the unit is not re-assembled properly, it will no longer function properly.
- Do not remove the retaining ring from the bag neck.

6.5 Manual cleaning and disinfection

- Solely manual reprocessing is not recommended due to the lower effectiveness and reproducibility. It should only be used if no automated procedure is available.
- Use fully demineralised water.
- Cleaning and disinfecting solution must be prepared before each manual cleaning cycle. Please follow the manufacturer's instructions for the solution.
- A solution of Sekusept Aktiv is recommended.
- Place the product in the cleaning and disinfecting solution. Pay attention to the manufacturer's recommended concentration and soaking time.
- Even when cleaning manually, make sure that all cavities are properly rinsed through.
- A small brush with soft bristles or a soft sponge can be used to help clean the product.
- After cleaning, rinse sufficiently with fully demineralised water to remove any residues.
- Avoid any water accumulation.
- Dry the product completely.
- Check visually for residues. Repeat the rinsing process if necessary.

6.6 Automatic cleaning and disinfection

- If using alkaline cleaning agents, neutralisation must be carried out first.
- Do not use any desiccants.

- It is recommended to use an automated/ mechanical cleaning and disinfection procedure.
- The product is suitable for machine processing/thermal disinfection ($\leq 93^{\circ}\text{C}$).
- It must be ensured that the selected washer-disinfector complies with ISO 15883.
- Follow the operating instructions of the respective device model.
- Use cleaning agents that are suitable for stainless steel and polysulfone.
- Recommended cleaning agents:
Cleaner: Sekumatic ProClean // Neutraliser: ECOLAB Sekumatic FNZ
- The operating cycle (cleaning, disinfection, rinsing, drying if necessary) must achieve the specified performance.
- Do not overload the basket and avoid shadowing.
- Load products with openings and recesses with open side facing downwards to avoid accumulation of cleaning water.

Instructions:

1. Place the products in the washer-disinfector/arrange around the spiral rack (ensure that all cavities are effectively flushed through). Observe the loading pattern of your washer-disinfector if necessary.
2. Start a cycle with the following parameters: (e.g. VARIO TD)
 - Pre-rinse with fully demineralised water 20°C // 1–4 min.
 - Cleaning with alkaline cleaner max. 55°C // 5–10 min.
 - Neutralisation 20°C // 2–6 min.
 - Intermediate rinsing 20°C // 2–3 min
 - Disinfection 93°C // 5 min.
 - Drying max. 100°C // 30 min
3. Inspect for any residual soil. Repeat the cycle if necessary.
4. Inspect the products and prepare for sterilization.

6.7 Drying

- After cleaning, the products must be dried thoroughly (preferably supported by compressed air).
- Do not use drying agents or cloths which could leave residues on the silicone.

6.8 Maintenance, control and inspection

- Inspect the products after disinfection and before each use for intactness, leakage and cleanliness.
- Check the product for cracks and material changes.
- If residual contamination is still visible, repeat the entire procedure until no more residues are visible.
- Reassemble the breathing bag as described in chapter 10.
- After reprocessing and before each use, a functional check must be carried out in accordance with the instructions in chapter 4 (at least item no 4.1).

6.9 Packaging

- After processing, package the product according to ISO 7153.
- Store sterilised products in sterile packaging in a closed cabinet, protected from dust, light, moisture and temperature fluctuations.

6.10 Sterilisation

- The product must be sterilised before each use.
- Before sterilisation, pack the products in a suitable container or in suitable sterilisation wrapping.
- The correct size of container or wrapping must be chosen.
- Use a sterilisation indicator for the wrapping and note the sterilisation and expiry date on the wrapping.
- The product is only suitable for steam sterilization in a fractionated vacuum process.
Maximum exposure time: 5 min.
Maximum temperature: 134°C
- Follow the instructions of the autoclave manufacturer carefully.
- After steam sterilisation, the product must be allowed to cool down at room temperature.
- The manufacturer accepts no liability for risks arising from improper sterilisation (contamination).

6.11 Additional Information

- Frequent sterilisation can accelerate the natural ageing of the product. In addition, the service life of sterilised products depends very much on carefully performed pre-treatment (cleaning and disinfection).
- When sterilising several products in one sterilisation cycle, the maximum load of the steriliser must not be exceeded.

7 Accessories and replacement parts

- The following three models are available.
The selection is according to the weight of the patient.
Observe the overview given in chapter 10.
- All models consist at least of patient valve, self-filling silicone breathing bag and inlet valve.
- Models with the suffix -OM, are supplied without a mask.
- Models with the suffix -OD, are supplied without a pressure limiting valve.
- A PEEP valve can be connected to all models in combination with an adapter.
- In case of compatibility questions please contact customer service.
- Only use replacement parts and accessories provided by Droh.
Otherwise the performance of the products may be limited.

Article no.

01-0M	Silicone mask round, size 0 (neonates)
01-1M	Silicone mask round, size 1 (infants)
01-2M	Silicone mask round, size 2 (toddlers)
01-3M	Transparent mask, silicone cushion, size 3 (children)
01-4M	Transparent mask, silicone cushion, size 4 (small adults)
01-5M	Transparent mask, silicone cushion, size 5 (adult)
01-05H	Patient valve expiration membrane, blue
01-06H	Patient valve lower part
01-07H	Patient valve lip membrane, blue
01-08H	Patient valve upper part with connection for pressure limiting valve
01-08H-OD	Patient valve upper part without connection for pressure limiting valve
01-09H	Silicone sealing ring blue for connection ring
01-10H	Connection ring for patient valve with sealing ring
01-11H	Adult bag (1700 ml), with connection ring
01-12H	Inlet valve Inlet membrane, blue
01-13H	Inlet valve inner part
01-14H	Inlet valve outer component with inlet socket and O ₂ connection nipple
01-15H	O ₂ -reservoir system for adults and children, complete with valve and bag (2700 ml)
01-19H	O ₂ reservoir bag (2700 ml) for adults and children
01-20H	Toddler bag (300 ml), with connection ring
01-21H	O ₂ reservoir bag(900 ml) for infants
01-22H	Patient valve complete (incl. pressure limiter) 01-22H-OD
01-22H-OD	Patient valve complete (excl. pressure limiter)
01-23H	Inlet valve complete
01-24H	O ₂ Reservoir valve complete
01-25H	Children's bag (470 ml), with connection ring
01-29H	Patient valve pressure limiting valve (60 mbar) single
01-31H	Patient valve pressure limiting valve (40 mbar) single
01-32H	Silicone hand strap for adult breathing bag
01-40H	Plastic case suitable for breathing bag and accessories
01-50H	O ₂ -Reservoir system for toddlers, complete with valve and bag (900 ml)

- For a complete list of compatible accessories please contact our Customer service or at **www.droh.de**.

8 Repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.
- If the cause of the error cannot be eliminated, the product must not be put into operation. After consulting the customer service, the product must be disposed of according to chapter 9 or returned to the sales distributor.
- In the event of a return for repair or replacement, the entire reprocessing process must be carried out beforehand in order to exclude any risk to the manufacturer's employees.
We reserve the right to reject soiled or contaminated products for safety reasons.

9 Disposal

- If any of the above-mentioned errors have occurred, the product must be disposed of after consulting customer service.
- The product must be disposed of at the latest **two years after first use or following 30 reprocessing cycles.**
- Safe disposal is carried out in accordance with current national regulations.

10 Technical Specifications

- The data recorded and listed as part of the specifications are mean values and may vary owing to changed conditions (e.g. temperature fluctuations/different levels of strength exerted by users when squeezing the breathing bag etc.).

		Infant breathing bags	Child breathing bags	Adult breathing bags
Patient weight		up to 8 kg	8 kg to 30 kg	over 30 kg
Bag	Weight	170 g	215 g	390 g
	Volume	300 ml	470 ml	1700 ml
Dimensions (L x W x H)		256 x 85 x 74 mm	255 x 91 x 91 mm	325 x 130 x 130 mm
Reservoir volume		900 ml	2700 ml	2700 ml
Pressure limiting system		40 cm H ₂ O (4,0 kPa)	40 cm H ₂ O (4,0 kPa)	60 cm H ₂ O (6,0 kPa)
Max. frequency / min		100	96	45
Max. tidal volume / ml		130	250	900
Inspiratory resistance			3.30 mbar at 50 l/m	
Expiratory resistance			2.20 mbar at 50 l/m	
Patient valve dead space			7 ml	
Patient connection		15 mm tapered socket part and 22 mm OD tapered plug part according to ISO 5356		
Expiration part		30 mm tapered plug part according to ISO 5356 on Diverter / Peep-Valve-Adapter		
Connection for oxygen		delivery nozzle 6 mm		
Connection for reservoir bag		6 mm OD/22 mm ID, conical connector		
Operating temperature		-20°C to +45°C		

Materials

Rigid plastic parts	Polysulfone
Soft plastic parts, valve membranes	Silicone rubber
Pressure limiting valve spring	Stainless steel
Reservoir bag	Med. polyvinylchloride (PVC)

■ **Inspiratory oxygen concentration**

Oxygen concentration depends on the configured oxygen flow, the ventilation volume and the ventilation frequency.

The following oxygen concentrations are achieved:

Breathing bag	Ventilation displacement x Ventilation frequency/min	Oxygen concentration in % with reservoir/(without reservoir) with an oxygen flow of litres per minute:					
		3 LPM	4 LPM	5 LPM	10 LPM	12 LPM	15 LPM
Infant breathing bag with 900 ml reservoir	20 x 40		97 (97)				
	40 x 30		97 (89)				
Child breathing bag with 1700 ml reservoir	100 x 30				97 (90)		
	250 x 20				97 (75)		
Adult breathing bag with 1700 ml reservoir	500 x 12	58 (37)		81 (52)	97 (73)	97 (84)	97 (89)
	500 x 24	39 (32)		52 (38)	84 (48)	97 (53)	97 (59)
	750 x 12	47 (33)		62 (41)	97 (56)	97 (62)	97 (69)
	750 x 24	34 (29)		42 (33)	65 (42)	74 (45)	86 (48)
	1000 x 12	41 (32)		52 (39)	84 (55)	94 (50)	97 (69)
	1000 x 24	30 (28)		38 (31)	53 (39)	59 (42)	66 (44)

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