Wolfram Droh GmbH

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REV GA/ 2021-04

Ningbo Huakun Medical Equipment Co., Ltd.
No. 208 Binjiang Road
Cixi High-tech Industrial Development Zone
Cixi City
315301 Zhejiang, P.R. China

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GB Contents

1	General	Page	9
2	Safety and Responsibilities	Page	9
3	Modification of the product	Page	10
4	Before using the product	Page	10
5	Using the product	Page	10
6	Product Maintenance / Reprocessing	Page	11
7	Accessories/Package content	Page	11
8	Replacement parts and repairs	Page	11
9	Disposal	Page	11
10	Explanation of symbols	Page	12

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1 General

- The product may only be used by trained professionals in cardiopulmonary reanimation and artificial ventilation.
- The anaesthetic masks are intended to be used in combination with a ventilator to ventilate a patient.
- The product is mainly used in hospitals and in emergency care.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- For traceability purposes, the product may only be stored in its original packaging. If this is not possible for the user, the user is responsible for correct labelling and traceability.
- If in doubt about the use and connection of the product, please contact customer service.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- Keep the product free of oil and grease, as this could present an explosion hazard. This includes the exposure to cleaning alcohols, hand cream and adhesives – in contact with highly compressed oxygen explosive reactions may occur.
- Only use replacement parts and accessories provided by Droh. Otherwise the performance of the products may be reduced.
- **Single use products may not be reused!**
They must be disposed of immediately after use in accordance with chapter 9.
The reuse of single use products represents a potential risk for the patient or the user. It can result in contamination and/or reduced functionality.
- Contamination and/or limited functionality of the product can lead to injury, illness or death of the patient.

3 Modification of the product

- Apart from the described maintenance instructions, the product must not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- Do not expose the product to direct sunlight or other light sources. Store and transport the product in a dry place in its original packaging at -40°C to $+70^{\circ}\text{C}$.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Verify that parts or components do not show any damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ Damaged products must not be used!
After consulting customer service, the product must be disposed of according to chapter 9.
- Please note that the product is supplied in a non-sterile condition.
- The hygienic preparation of the product lies within the clinical responsibility of the user! The plastic product can be disinfected with a suitable cleaning agent (e.g. alcohol) before use. Sterilization is not permitted!
- Attach the mask to a suitable ventilator.
- The product is now ready for use.

5 Using the product

- The operating temperature should be between -40°C and $+50^{\circ}\text{C}$.
- Place the mask over the patient's face, so that no air can escape from the mask.
- Depending on the model a cannula (Luer connector) can be used to adjust the inflation of the cushioned rim.
- Choose the correct type and size of mask according to the anatomical features of the patient.
If applicable observe national legislation.

6 Product maintenance / Reprocessing

- Single use products may not be reused!
- Please discard the product after use in accordance with the instructions in chapter 9.

7 Accessories / Package content

- In case of compatibility questions please contact customer service.

8 Replacement parts and repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.
- If the cause of the error cannot be eliminated, the product must not be put into operation. After consulting the customer service, the product must be disposed of according to chapter 9 or returned to the sales distributor.

9 Disposal

- The product must be disposed of after use.
- Please note the expiry date if applicable.
The product must not be used after the expiry date!
- Safe disposal is carried out in accordance with current national regulations.

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